

Anti-TM7SF3 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A15928-1-Fluoro550

About TM7SF3

Involved in cellular response to unfolded protein; negative regulation of programmed cell death; and positive regulation of insulin secretion. Located in plasma membrane.

Overview

Product Name	Anti-TM7SF3 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NS93

Technical Details

Immunogen	E.coli-derived human TM7SF3 recombinant protein (Position: R102-R544).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TM7SF3 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.