

Anti-LEKR1 Antibody Picoband® Cy3 Conjugated

Catalog Number: A15972-1-Cy3

About LEKR1

LEKR1 (Leucine, Glutamate And Lysine Rich 1, also known as FLJ16641 or FLJ37161, LOC389170) is located on human chromosome 3q25.31. The LEKR1 protein is biasedly expressed in testis, thyroid and other tissues. Structurally, the encoded protein is reported to be 45158 Da in mass. LEKR1 protein is expected to localize in various compartments (mitochondrion, nucleus). LEKR1 is associated with Obsessive-Compulsive Personality Disorder.

Overview

Product Name	Anti-LEKR1 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	J3KP02

Technical Details

Immunogen	E.coli-derived human LEKR1 recombinant protein (Position: K590-R688).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LEKR1 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.