

Anti-ACTL9 Antibody

Catalog Number: A16009-1

About ACTL9

ACTL9 (Actin Like 9), also named MGC33407, encodes a testis-specific actin-like protein that plays a role in fusion of proacrosomal vesicles and perinuclear theca formation (PMID: 33626338). ACTL9 protein contains 416 amino acids, which mainly localizes in the equatorial segment of human sperm head (PubMed: 33626338). Diseases associated with ACTL9 include spermatogenic failure and non-syndromic male infertility.

Overview

Product Name	Anti-ACTL9 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ACTL9 Antibody catalog # A16009-1. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8TC94

Technical Details

Immunogen	E.coli-derived human ACTL9 recombinant protein (Position: A3-S169).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ACTL9 Antibody

For Research Use Only. Not for use in diagnostic procedures.