

Anti-Brain protein I3 BRI3 Antibody

Catalog Number: A16038

About BRI3

Belongs to the glycosyl hydrolase 35 family.

Ota T., Nat. Genet. 36:40-45(2004).

Taylor T.D., Nature 440:497-500(2006).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Overview

Product Name	Anti-Brain protein I3 BRI3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Brain protein I3 BRI3 Antibody catalog # A16038. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95415

Technical Details

Immunogen	Synthesized peptide derived from internal of human GLB1L3.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.
Some PubMed article(s) citing the expression level of this target are as follows:
Boster Bio's internal QC testing used:
WB: 1:500-1:1000

Anti-Brain protein I3 BRI3 Antibody (A16038) Images




Figure 1. Western blotting validation for Anti-Brain protein I3 BRI3 Antibody A16038

Western blot (WB) analysis of BRI3 polyclonal antibody at 1:500 dilution

Lane1:Hela whole cell lysate

Lane2:Mouse brain tissue lysate

Lane3:Rat colon tissue lysate

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Brain protein I3 BRI3 Antibody