

Anti-Mob3A/B Antibody

Catalog Number: A16079


Overview

Product Name	Anti-Mob3A/B Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Mob3A/B Antibody catalog # A16079. Tested in IHC, WB applications. This antibody reacts with Human, Mouse.
Application	IHC, WB
Clonality	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96BX8


Technical Details

Immunogen	Synthesized peptide derived from the N-terminal region of human Mob3A/B. at AA range: 1-80
Isotype	IgG
Form	Liquid
Purification	Immunogen affinity purified
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB, 1:500-1:2000 IHC-P, 1:100-1:300</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-Mob3A/B Antibody (A16079) Images



Western Blot (WB) analysis of 3T3 using Mob3A/B antibody.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Skin, antibody was diluted at 1:200.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Mob3A/B Antibody