

Anti-TBC1D3A/B/C TBC1D3C Antibody

Catalog Number: A16395

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-TBC1D3A/B/C TBC1D3C Antibody
Reactive Species	Human
Description	Boster Bio Anti-TBC1D3A/B/C TBC1D3C Antibody catalog # A16395. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6IPX1


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human TBC1D3. AA range:481-530
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Suggested Dilutions

WB 1:500-1:2000
ELISA 1:20000
IF 1:50-200

Anti-TBC1D3A/B/C TBC1D3C Antibody (A16395) Images



Western blotting validation for Anti-TBC1D3A/B/C TBC1D3C Antibody A16395 Western Blot (WB) analysis of specific cells using TBC1D3A/B/C polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Anti-TBC1D3A/B/C TBC1D3C Antibody

For Research Use Only. Not for use in diagnostic procedures.