

Anti-ZNF90 Antibody (Center)

Catalog Number: A16460

About ZNF90

ZNF90 may be involved in transcriptional regulation.


Overview

Product Name	Anti-ZNF90 Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-ZNF90 Antibody (Center) (Catalog # A16460). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q03938

Technical Details

Immunogen	This ZNF90 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 84-112 amino acids from the Central region of human ZNF90.
Predicted Reactive Species	Rat
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-ZNF90 Antibody (Center) (A16460) Images



ZNF90 Antibody (Center) western blot analysis in CEM cell line lysates (35ug/lane). This demonstrates the ZNF90 antibody detected the ZNF90 protein (arrow).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZNF90 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.