

## Anti-Olfactory receptor 2W1 OR2W1 Antibody

Catalog Number: A16478

### Overview

Product Name	Anti-Olfactory receptor 2W1 OR2W1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Olfactory receptor 2W1 OR2W1 Antibody catalog # A16478. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y3N9

### Technical Details

Immunogen	Synthesized peptide derived from part region of human protein AA range: 122-172
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	OR2W1 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Olfactory receptor 2W1 OR2W1 Antibody

For Research Use Only. Not for use in diagnostic procedures.