

## Anti-Olfactory receptor 10H1 Antibody

Catalog Number: A16480

### Overview

Product Name	Anti-Olfactory receptor 10H1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Olfactory receptor 10H1 Antibody catalog # A16480. Tested in WB, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9Y4A9


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human OR10H1. AA range:236-285
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:20000

## Anti-Olfactory receptor 10H1 Antibody (A16480) Images



Western Blot analysis of various cells using Olfactory receptor 10H1 Polyclonal Antibody



Western blot analysis of lysates from HeLa and MCF-7 cells, using OR10H1 Antibody. The lane on the right is blocked with the synthesized peptide.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Olfactory receptor 10H1 Antibody

For Research Use Only. Not for use in diagnostic procedures.