

Anti-CEP170/CEP170P1 Antibody

Catalog Number: A16502

About CEP170P1

This locus appears to be a transcribed pseudogene similar to centrosomal protein 170kDa (CEP170). An approximately 50 kb region upstream of this locus also is homologous to CEP170, but is not transcribed.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-CEP170/CEP170P1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CEP170/CEP170P1 Antibody catalog # A16502. Tested in WB, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96L14

Technical Details

Immunogen	E.coli-derived human CEP170/CEP170P1 recombinant protein (Position: T3-T290).
Form	Liquid
Concentration	500 ug/ml

Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400

[Anti-CEP170/CEP170P1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.