

Anti-Olfactory receptor 10T2 Antibody

Catalog Number: A16511

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Olfactory receptor 10T2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Olfactory receptor 10T2 Antibody catalog # A16511. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NGX3


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human OR10T2. AA range:231-280
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified


Suggested Dilutions

WB 1:500-1:2000
ELISA 1:20000

Anti-Olfactory receptor 10T2 Antibody (A16511) Images



Western Blot analysis of 293T-UV cells using Olfactory receptor 10T2 Polyclonal Antibody diluted at 1:500



Western blot analysis of lysates from HUVEC and COLO cells, using OR10T2 Antibody. The lane on the right is blocked with the synthesized peptide.

[Anti-Olfactory receptor 10T2 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.