

Anti-IGLON5 Antibody

Catalog Number: A16589


Overview

Product Name	Anti-IGLON5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IGLON5 Antibody catalog # A16589. Tested in WB applications. This antibody reacts with Human,Rat,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	A6NGN9

Technical Details

Immunogen	Synthetic peptide, corresponding to Human IGLON5.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:2000

Anti-IGLON5 Antibody (A16589) Images



Western blot (WB) analysis of IGLON5 polyclonal antibody at 1:1000 dilution Lane1:BV2 whole cell lysate(40ug) Lane2:A2780 whole cell lysate(40ug) Lane3:EC9706 whole cell lysate(40ug) Lane4:C6 whole cell lysate(20ug) Lane5:U-87MG whole cell lysate(30ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IGLON5 Antibody

For Research Use Only. Not for use in diagnostic procedures.