

Anti-LCN12 Antibody

Catalog Number: A16605-1

About LCN12

LCN12 (Lipocalin 12) is an epididymis-specific protein that might play a significant physiological role in male reproduction (PMID: 20692383).

Overview

Product Name	Anti-LCN12 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-LCN12 Antibody catalog # A16605-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6JVE5

Technical Details

Immunogen	E.coli-derived human LCN12 recombinant protein (Position: F53-D174).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LCN12 Antibody

For Research Use Only. Not for use in diagnostic procedures.