

Anti-hnRNP CL1 Antibody

Catalog Number: A16614

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-hnRNP CL1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-hnRNP CL1 Antibody catalog # A16614. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60812

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human HNRCL. AA range:31-80
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.
Suggested Dilutions	WB 1:500-1:2000

	ELISA 1:20000
--	---------------

Anti-hnRNP CL1 Antibody (A16614) Images



Western blot analysis of HNRCL Antibody. The lane on the right is blocked with the HNRCL peptide.

Anti-hnRNP CL1 Antibody

For Research Use Only. Not for use in diagnostic procedures.