

Anti-RB40A RAB40A Antibody

Catalog Number: A16646

About RAB40A

2 isoforms of the human protein are produced by alternative splicing.

Bechtel S., BMC Genomics 8:399-399(2007).

Nusbaum C., Nature 437:551-555(2005).

Ota T., Nat. Genet. 36:40-45(2004).

Overview

Product Name	Anti-RB40A RAB40A Antibody
Reactive Species	Human
Description	Boster Bio Anti-RB40A RAB40A Antibody catalog # A16646. Tested in ELISA, WB applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WXH6

Technical Details

Immunogen	Synthesized peptide derived from human RB40A protein.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	RB40A antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.
Some PubMed article(s) citing the expression level of this target are as follows:
Boster Bio's internal QC testing used:
WB 1:500-2000 ELISA 1:5000-20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RB40A RAB40A Antibody