

Anti-SPRR2F Antibody

Catalog Number: A16672

Overview

Product Name	Anti-SPRR2F Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-SPRR2F Antibody catalog # A16672. Tested in IHC, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96RM1

Technical Details

Immunogen	E.coli-derived human SPRR2F recombinant protein (Position: P19-Q51).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPRR2F Antibody

For Research Use Only. Not for use in diagnostic procedures.