

Anti-AN1-type zinc finger protein 3 ZFAND3 Antibody

Catalog Number: A16734


Overview

Product Name	Anti-AN1-type zinc finger protein 3 ZFAND3 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-AN1-type zinc finger protein 3 ZFAND3 Antibody catalog # A16734. Tested in IHC applications. This antibody reacts with Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H8U3


Technical Details

Immunogen	Recombinant fusion protein of human ZFAND3(NP_068762.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-AN1-type zinc finger protein 3 ZFAND3 Antibody (A16734) Images



Immunohistochemistry of paraffin-embedded rat kidney using ZFAND3 antibody at dilution of 1:100 .Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded mouse kidney using ZFAND3 antibody at dilution of 1:100 .Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AN1-type zinc finger protein 3 ZFAND3 Antibody

For Research Use Only. Not for use in diagnostic procedures.