

Anti-Olfactory receptor 5111 Antibody

Catalog Number: A16884

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Olfactory receptor 5111 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Olfactory receptor 5111 Antibody catalog # A16884. Tested in ELISA, IF, WB applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H343


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human OR5111. AA range:265-314
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Suggested Dilutions

WB 1:500-1:2000
IF 1:200-1:1000
ELISA 1:10000

Anti-Olfactory receptor 5111 Antibody (A16884) Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night



Immunofluorescence analysis of COS7 cells, using OR5111 Antibody. The picture on the right is blocked with the synthesized peptide.

Anti-Olfactory receptor 5111 Antibody

For Research Use Only. Not for use in diagnostic procedures.