

Anti-GLP-1/ZGLP1 Antibody

Catalog Number: A17057-1

About ZGLP1

ZGLP1, also known as GLP-1, is a nuclear protein containing two zinc fingers, one of which is highly homologous to the conserved GATA zinc finger DNA binding motif. It is a transcriptional repressor that plays a central role in somatic cells of the gonad and is required for germ cell development.

Overview

Product Name	Anti-GLP-1/ZGLP1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-GLP-1/ZGLP1 Antibody catalog # A17057-1. Tested in WB, IHC, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	POC6A0

Technical Details

Immunogen	Peptide
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GLP-1/ZGLP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.