

Anti-ARPM1/ACTRT3 Antibody

Catalog Number: A17256-2

About ACTRT3

Predicted to be located in cytoplasm; cytoskeleton; and male germ cell nucleus. Predicted to be active in actin cytoskeleton.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-ARPM1/ACTRT3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ARPM1/ACTRT3 Antibody catalog # A17256-2. Tested in WB, IHC, ICC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9BYD9

Technical Details

Immunogen	E.coli-derived human ARPM1/ACTRT3 recombinant protein (Position: K18-P230).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.

Suggested Dilutions

Western blot, 1:500-2000
Immunohistochemistry, 1:50-400
Immunocytochemistry/Immunofluorescence, 1:50-400

Anti-ARPM1/ACTRT3 Antibody

For Research Use Only. Not for use in diagnostic procedures.