

## Anti-Calpain-12 CAPN12 Antibody

Catalog Number: A17322-1


### Overview

Product Name	Anti-Calpain-12 CAPN12 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Calpain-12 CAPN12 Antibody catalog # A17322-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6ZSI9

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 222-266 of Human CAPN12.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 96% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-Calpain-12 CAPN12 Antibody (A17322-1) Images



Western blot (WB) analysis of CAPN12 polyclonal antibody at 1:500 dilution Lane1:A549 whole cell lysate(40ug) Lane2:Jurkat whole cell lysate(40ug) Lane3:THP-1 whole cell lysate(40ug) Lane4:The lung tissue lysate of Mouse(40ug) Lane5:The lung tissue lysate of Fetal Rat(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Calpain-12 CAPN12 Antibody

For Research Use Only. Not for use in diagnostic procedures.