

Anti-Actin-related protein 10 ACTR10 Antibody

Catalog Number: A17479

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Actin-related protein 10 ACTR10 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Actin-related protein 10 ACTR10 Antibody catalog # A17479. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NZ32

Technical Details

Immunogen	Recombinant fusion protein of human ACTR10(NP_060947.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-Actin-related protein 10 ACTR10 Antibody (A17479) Images



Anti-Actin-related protein 10 ACTR10 Antibody

For Research Use Only. Not for use in diagnostic procedures.