

Anti-SNX32 Antibody Picoband® HRP Conjugated

Catalog Number: A17611-1-HRP

About SNX32

Predicted to enable phosphatidylinositol binding activity. Predicted to be involved in retrograde transport, endosome to Golgi. Predicted to be located in cytoplasmic vesicle and cytosol. Predicted to be active in endosome.

Overview

Product Name	Anti-SNX32 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q86XE0

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human SNX32.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-SNX32 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.