

Anti-OR2T3/2T34 Antibody

Catalog Number: A17737

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-OR2T3/2T34 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-OR2T3/2T34 Antibody (Catalog# A17737). Tested in WB, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, ICC, WB
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NH03

Technical Details


Immunogen	A peptide derived from human OR2T3/2T34. Immunogen sequence location: 269 - 318
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	WB 1:500 - 1:2000.

	ICC/IF 1:200 - 1:1000. ELISA: 1:10000. Not yet tested in other applications.
--	--


Anti-OR2T3/2T34 Antibody (A17737) Images



Western blot analysis of the lysates from 293 cells using OR2T3/34 antibody.



Immunofluorescence analysis of A549 cells, using OR2T3/2T34 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from MCF-7, HUVEC, HeLa, and HepG2 cells, using OR2T3/2T34 Antibody. The lane on the right is blocked with the synthesized peptide.

Anti-OR2T3/2T34 Antibody

For Research Use Only. Not for use in diagnostic procedures.