

## Anti-Tubulin alpha-3C/D/E Antibody

Catalog Number: A17842


### Overview

Product Name	Anti-Tubulin alpha-3C/D/E Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Tubulin alpha-3C/D/E Antibody catalog # A17842. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13748


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human TUBA3C/E. AA range:201-250
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:5000

## Anti-Tubulin alpha-3C/D/E Antibody (A17842) Images



Western Blot analysis of various cells using Tubulin alpha-3C/D/E Polyclonal Antibody diluted at 1:2000. Secondary antibody was diluted at 1:20000



Western Blot analysis of HEPG2 HELA cells using Antibody diluted at 2000. Secondary antibody was diluted at 1:20000

## 2 Publications Citing This Product

1. PubMed ID: 10.1007/s12192-008-0086-5, Prohibitin protects against oxidative stress-induced cell injury in cultured neonatal cardiomyocyte
2. PubMed ID: 23326455, Meng R, Zhu D, Bi Y, Yang D, Wang Y. Plos One. 2013;8(1):E53557. Doi: 10.1371/Journal.Pone.0053557. Epub 2013 Jan 10. Erythropoietin Inhibits Gluconeogenesis And Inflammation In The Liver And Improves Glucose Intolerance In High-Fat Diet-Fed Mice.

Visit [bosterbio.com/anti-tubulin-alpha-3c-d-e-antibody-a17842-boster.html](http://bosterbio.com/anti-tubulin-alpha-3c-d-e-antibody-a17842-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tubulin alpha-3C/D/E Antibody

For Research Use Only. Not for use in diagnostic procedures.