

Anti-FAM127/RTL8A Antibody

Catalog Number: A17938

About RTL8A

The FAM127 family currently has three members, namely FAM127A (also called RTL8C and CXX1), FAM127B (also called CXX1B and RTL8A), and FAM127C (also called RTL8B and CXX1C). They are located on chromosomes and their functions are unknown. All three are composed of 113 amino acids and have extremely high sequence similarity. The immunogen of this antibody has a sequence similarity of more than 95% with members of the FAM127 family. Theoretically, the antibody can recognize all members of the FAM127 family.

Overview

Product Name	Anti-FAM127/RTL8A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-FAM127/RTL8A Antibody catalog # A17938. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9BWD3

Technical Details

Immunogen	E.coli-derived human FAM127/RTL8A recombinant protein (Position: G3-E111).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300

ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FAM127/RTL8A Antibody

For Research Use Only. Not for use in diagnostic procedures.