

## Anti-LRRC71 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A17988-1-Fluoro488

### About LRRC71

The function of this protein remains unknown.

### Overview

Product Name	Anti-LRRC71 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8N4P6

### Technical Details

Immunogen	E.coli-derived human C1orf92/LRRC71 recombinant protein (Position: K35-D543). Human C1orf92/LRRC71 shares 81.6% amino acid (aa) sequence identity with mouse C1orf92/LRRC71.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LRRC71 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.