

## Anti-LRRC71 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A17988-1-Fluoro550

### About LRRC71

The function of this protein remains unknown.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-LRRC71 Antibody Picoband® Fluoro550 Conjugated  |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | Q8N4P6   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human C1orf92/LRRC71 recombinant protein (Position: K35-D543). Human C1orf92/LRRC71 shares 81.6% amino acid (aa) sequence identity with mouse C1orf92/LRRC71. |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | Fluoro550<br>Excitation Wavelength: 562 nm<br>Emission Wavelength: 576 nm  |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.   |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LRRC71 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.