

Anti-LRRC71 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A17988-1-Fluoro647

About LRRC71

The function of this protein remains unknown.

Overview

| | |
|----------------------|--|
| Product Name | Anti-LRRC71 Antibody Picoband® Fluoro647 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q8N4P6 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human C1orf92/LRRC71 recombinant protein (Position: K35-D543). Human C1orf92/LRRC71 shares 81.6% amino acid (aa) sequence identity with mouse C1orf92/LRRC71. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LRRC71 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.