

Anti-MFSD13A Antibody Picoband® APC Conjugated

Catalog Number: A17995-1-APC

About MFSD13A

MFSD13A (Major Facilitator Superfamily Domain Containing 13A) is located on human chromosome 10q24.32. The MFSD13A protein is ubiquitously expressed in brain, lung and other tissues. Structurally, the encoded protein is reported to be 27160 Da in mass. The MFSD13A (also known as TMEM180; C10orf77; bA18l14.8) gene is conserved in chimpanzee, Rhesus monkey, dog, cow, mouse, chicken, frog, etc. 234 organisms have orthologs with human gene MFSD13A. The MFSD13A gene might be genetic modifiers of pancreatic cancer tumorigenesis.

Overview

Product Name	Anti-MFSD13A Antibody Picoband® APC Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q14CX5

Technical Details

Immunogen	E.coli-derived human MFSD13A recombinant protein (Position: N66-V517).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MFSD13A Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.