

## Anti-ZMYM1 Antibody Picoband® APC Conjugated

Catalog Number: A18075-1-APC

### About ZMYM1

Zinc finger MYM-type containing 1 is a protein that in humans is encoded by the ZMYM1 gene. Predicted to enable protein dimerization activity and zinc ion binding activity. Predicted to be located in nucleus.

### Overview

Product Name	Anti-ZMYM1 Antibody Picoband® APC Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5SVZ6

### Technical Details

Immunogen	E.coli-derived human ZMYM1 recombinant protein (Position: R454-E1124).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ZMYM1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.