

Anti-ZMYM1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A18075-1-Fluoro488

About ZMYM1

Zinc finger MYM-type containing 1 is a protein that in humans is encoded by the ZMYM1 gene. Predicted to enable protein dimerization activity and zinc ion binding activity. Predicted to be located in nucleus.

Overview

Product Name	Anti-ZMYM1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5SVZ6

Technical Details

Immunogen	E.coli-derived human ZMYM1 recombinant protein (Position: R454-E1124).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZMYM1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.