

Anti-TTC36 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A18139-Fluoro594

About TTC36

The protein encoded by this gene has three tetratricopeptide repeats and is a chaperone for heat shock protein 70. The encoded protein may function as a tumor suppressor in hepatocellular carcinoma (HCC) since it promotes apoptosis but is downregulated in HCC.

Overview

Product Name	Anti-TTC36 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A6NLP5

Technical Details

Immunogen	E.coli-derived human TTC36 recombinant protein (Position: M1-D187). Human TTC36 shares 89.2% and 88.6% amino acid (aa) sequence identity with mouse and rat TTC36, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TTC36 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.