

Anti-SPINK13 Antibody Picoband® Cy3 Conjugated

Catalog Number: A18284-Cy3

About SPINK13

SPINK13 cDNA/gene is a gene with protein product which located on 5q32. The SPINK13 gene is conserved in chimpanzee, Rhesus monkey, mouse, and rat. 65 organisms have orthologs with human gene SPINK13.

Overview

| | |
|----------------------|--|
| Product Name | Anti-SPINK13 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q1W4C9 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human SPINK13 recombinant protein (Position: I25-D94). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPINK13 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.