

Anti-SPINK13 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A18284-Fluoro647

About SPINK13

SPINK13 cDNA/gene is a gene with protein product which located on 5q32. The SPINK13 gene is conserved in chimpanzee, Rhesus monkey, mouse, and rat. 65 organisms have orthologs with human gene SPINK13.

Overview

Product Name	Anti-SPINK13 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q1W4C9

Technical Details

Immunogen	E.coli-derived human SPINK13 recombinant protein (Position: I25-D94).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPINK13 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.