

Anti-SPINK13 Antibody Picoband® HRP Conjugated

Catalog Number: A18284-HRP

About SPINK13

SPINK13 cDNA/gene is a gene with protein product which located on 5q32. The SPINK13 gene is conserved in chimpanzee, Rhesus monkey, mouse, and rat. 65 organisms have orthologs with human gene SPINK13.

Overview

Product Name	Anti-SPINK13 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q1W4C9

Technical Details

Immunogen	E.coli-derived human SPINK13 recombinant protein (Position: I25-D94).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SPINK13 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.