

Anti-C16orf71 Antibody (Center)

Catalog Number: A18346

About DNAAF8

The exact function of C16orf71 remains unknown.

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-C16orf71 Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-C16orf71 Antibody (Center) (Catalog # A18346). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8IYS4

Technical Details

Immunogen	This C16orf71 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 320-348 amino acids from the Central region of human C16orf71.
Predicted Reactive Species	Bovine
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-C16orf71 Antibody (Center) (A18346) Images



C16orf71 Antibody (Center) western blot analysis in HL-60, MCF-7 and Raji cell line lysates (35ug/lane). This demonstrates the C16orf71 antibody detected the C16orf71 protein (arrow).

Anti-C16orf71 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.