

## Anti-KLHL28 Antibody

Catalog Number: A18441-1

### About KLHL28

KLHL28 also named as BTBD5 or Kelch like protein 28 is a 571 amino acid protein, which contains one BTB domain and six kelch repeats. KLHL28 localizes in the cytoplasm and belongs to BTB/kelch-associated proteins family. KLHL28 may play a role in organizing the actin cytoskeleton of the brain cells. KLHL28 is expressed as two isoforms produced by alternative splicing.

### Overview

Product Name	Anti-KLHL28 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-KLHL28 Antibody catalog # A18441-1. Tested in WB, IHC, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NXS3

### Technical Details

Immunogen	E.coli-derived human KLHL28 recombinant protein (Position: V224-L568).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-KLHL28 Antibody

For Research Use Only. Not for use in diagnostic procedures.