

Anti-LMNTD2 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A18456-Fluoro594

About LMNTD2

Predicted to be a structural constituent of chromatin. Predicted to act upstream of or within positive regulation of mRNA splicing, via spliceosome. Predicted to be active in lamin filament.

Overview

Product Name	Anti-LMNTD2 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8IXW0

Technical Details

Immunogen	E.coli-derived human LMNTD2 recombinant protein (Position: E155-A634). Human LMNTD2 shares 48.5% amino acid (aa) sequence identity with mouse LMNTD2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-LMNTD2 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.