

## Anti-OVOS1 Antibody (Center)

Catalog Number: A18488

### About OVOS1

Is able to inhibit all four classes of proteinases by a unique 'trapping' mechanism (By similarity).


### Overview

Product Name	Anti-OVOS1 Antibody (Center)
Reactive Species	Human
Description	Boster Bio Anti-OVOS1 Antibody (Center) (Catalog # A18488). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6IE37

### Technical Details

Immunogen	This OVOS1 antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 561-595 amino acids from the Central region of human OVOS1.
Predicted Reactive Species	Bovine
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

## Anti-OVOS1 Antibody (Center) (A18488) Images



Western blot analysis of lysate from human brain tissue lysate, using OVOS1 Antibody (Center). A18488 was diluted at 1:1000. A goat anti-rabbit IgG H&L (HRP) at 1:5000 dilution was used as the secondary antibody. Lysate at 35ug.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-OVOS1 Antibody (Center)

For Research Use Only. Not for use in diagnostic procedures.