

Anti-SLC35E4 Antibody Picoband® Biotin Conjugated

Catalog Number: A18522-Biotin

About SLC35E4

Predicted to enable antiporter activity. Predicted to be involved in transmembrane transport. Predicted to be integral component of membrane. Predicted to be active in Golgi apparatus.

Overview

Product Name	Anti-SLC35E4 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6ICL7

Technical Details

Immunogen	E.coli-derived human SLC35E4 recombinant protein (Position: M1-L350).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC35E4 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.