

Anti-SLC35E4 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A18522-Fluoro488

About SLC35E4

Predicted to enable antiporter activity. Predicted to be involved in transmembrane transport. Predicted to be integral component of membrane. Predicted to be active in Golgi apparatus.

Overview

Product Name	Anti-SLC35E4 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6ICL7

Technical Details

Immunogen	E.coli-derived human SLC35E4 recombinant protein (Position: M1-L350).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC35E4 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.