

## Anti-TOMM20L Antibody Picoband® PE Conjugated

Catalog Number: A18568-PE

### About TOMM20L

Predicted to enable mitochondrion targeting sequence binding activity. Predicted to contribute to protein transmembrane transporter activity. Predicted to be involved in protein import into mitochondrial matrix and tRNA import into mitochondrion. Predicted to be integral component of membrane. Predicted to be part of mitochondrial outer membrane translocase complex. Predicted to be integral component of mitochondrial outer membrane.

### Overview

Product Name	Anti-TOMM20L Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6UXN7

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human TOMM20L, which shares 70% amino acid (aa) sequence identity with mouse TOMM20L.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TOMM20L Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.