

Anti-RAB44 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A18664-Fluoro550

About RAB44

RAB44, member RAS oncogene family is a protein that in humans is encoded by the RAB44 gene. Predicted to enable GTP binding activity; GTPase activity; and calcium ion binding activity. Predicted to be located in azurophil granule membrane; plasma membrane; and specific granule membrane.

Overview

| | |
|----------------------|--|
| Product Name | Anti-RAB44 Antibody Picoband® Fluoro550 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q7Z6P3 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human RAB44 recombinant protein (Position: Q175-R995). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB44 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.