

## Anti-RAB44 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A18664-Fluoro594

### About RAB44

RAB44, member RAS oncogene family is a protein that in humans is encoded by the RAB44 gene. Predicted to enable GTP binding activity; GTPase activity; and calcium ion binding activity. Predicted to be located in azurophil granule membrane; plasma membrane; and specific granule membrane.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-RAB44 Antibody Picoband® Fluoro594 Conjugated   |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | Q7Z6P3   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E.coli-derived human RAB44 recombinant protein (Position: Q175-R995).     |
| Cross Reactivity    | No cross-reactivity with other proteins.                                  |
| Isotype             | Rabbit IgG  |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | Fluoro594<br>Excitation Wavelength: 593 nm<br>Emission Wavelength: 618 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.      |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB44 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.