

Anti-RAB44 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A18664-Fluoro647

About RAB44

RAB44, member RAS oncogene family is a protein that in humans is encoded by the RAB44 gene. Predicted to enable GTP binding activity; GTPase activity; and calcium ion binding activity. Predicted to be located in azurophil granule membrane; plasma membrane; and specific granule membrane.

Overview

Product Name	Anti-RAB44 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q7Z6P3

Technical Details

Immunogen	E.coli-derived human RAB44 recombinant protein (Position: Q175-R995).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB44 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.