

Anti-LRRC75B Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A18683-1-Fluoro550

About LRRC75B

FAM211B (family with sequence similarity 211, member B) is a leucine-rich repeat (LRR) protein present in viruses as well as eukaryotes. It is involved in diverse functions and is involved in protein-protein interactions. The LRR domains fold into a horseshoe (or arc) shape with a parallel beta-sheet on the concave face.

Overview

Product Name	Anti-LRRC75B Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q2VPJ9

Technical Details

Immunogen	E.coli-derived human LRRC75B recombinant protein (Position: Q59-E283).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LRRC75B Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.