

Anti-Olfactory receptor 10D4 OR10D4P Antibody

Catalog Number: A18752


Overview

Product Name	Anti-Olfactory receptor 10D4 OR10D4P Antibody
Reactive Species	Human
Description	Boster Bio Anti-Olfactory receptor 10D4 OR10D4P Antibody catalog # A18752. Tested in WB, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NGN7

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human OR10D4. AA range:161-210
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Anti-Olfactory receptor 10D4 OR10D4P Antibody (A18752) Images



Western Blot (WB) analysis of specific cells using Olfactory receptor 10D4 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Olfactory receptor 10D4 OR10D4P Antibody

For Research Use Only. Not for use in diagnostic procedures.