

Anti-SLCO1B7 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A18877-1-Fluoro550

About SLCO1B7

Predicted to enable bile acid transmembrane transporter activity and sodium-independent organic anion transmembrane transporter activity. Predicted to be involved in bile acid and bile salt transport and sodium-independent organic anion transport. Predicted to be located in plasma membrane. Predicted to be active in basolateral plasma membrane.

Overview

Product Name	Anti-SLCO1B7 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	G3V0H7

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human SLCO1B7.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-SLCO1B7 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.