

Anti-C2orf74 Antibody (C-term)

Catalog Number: A19009

About C2orf74

The function of this protein is unknown.


Overview

Product Name	Anti-C2orf74 Antibody (C-term)
Reactive Species	Human
Description	Boster Bio Anti-C2orf74 Antibody (C-term) (Catalog # A19009). Tested in WB application(s). This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide.
Storage Instructions	Maintain refrigerated at 2-8°C for up to 2 weeks. For long-term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Host	Rabbit
Uniprot ID	A8MZ97

Technical Details

Immunogen	This C2orf74 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 164-194 amino acids from the C-terminal region of human C2orf74.
Predicted Reactive Species	Bovine, Chicken, Mouse, Pig, Rat, Zebrafish
Isotype	Rabbit IgG
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Suggested Dilutions	WB: 1:1000

Anti-C2orf74 Antibody (C-term) (A19009) Images



C2orf74 Antibody (C-term) western blot analysis in human placenta tissue lysates (35ug/lane). This demonstrates the C2orf74 antibody detected the C2orf74 protein (arrow).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C2orf74 Antibody (C-term)

For Research Use Only. Not for use in diagnostic procedures.